


Wireless Multi-Channel Central Monitoring Unit



Monitor multiple residents from one central location

Eliminate In-Room Alarm Noise

The Smart Caregiver 433-CMU Central Monitoring Unit allows you to monitor multiple wireless components from one central location, such as a Nurses' Station. This eliminates in-room alarm noise & significantly reduces the tripping hazards associated with corded units.

Know Who Needs Help

The Central Monitor alerts audibly and the display visually indicates which resident needs assistance. At a glance, this provides caregivers with the information they need to protect their residents from falls.

Local Pager Notification

The CMU communicates with any number of 433-Pager Units, and upon activation, the pager audibly alerts and LCD screen displays which component activated the alert. This enables caregivers to attend to their duties yet be alerted in real time that a resident needs assistance. 433-PRB has a reset button that resets the 433-CMU-40/60 when alarming.

Key Features:

- Alarms audibly & visually at a central location and portable caregiver pager(s)
- Compatible with all Smart wireless components
- 40 channel capacity Range up to 300 feet
- Systems may be pre-programmed prior to delivery or easily re-programmed using the provided instructions
- AC adapter included (AC-04)
- Batteries not included (4-C)

Prevent Falls & Wandering

Easy-To-Install Call Button & Light System



**Arrives pre-programmed
and ready to use!**

Easy Installation

The 4017-SYS Emergency Call Light System includes a wireless call button with pull cord, reset button, and call alert with light. It is ideal for bathroom and hallway notification systems. All components are easily mounted on most surfaces with no hard wiring required.

4017-SYS complete system includes:

- 4017-CB: Wall Mount Call Button with Pull Cord
- 4017-RB: Wall Mount Reset Button
- 4017-ECA: Emergency Call Alert

Key Features:

Call Light

- Adjustable volume
- Audible & visible alert
- Works with up to 10 components
- No hard wiring
- Ideal Hallway Alert
- Batteries not included (3-C)
- Optional AC adapter available (AC-04)

Call Buttons

- Up to 300 feet
- Water Resistant
- Can be programmed to the TL-4015 Central Monitor
- No hard wiring
- Ideal for bathroom use

4017-SYS

Call Us Today! (888) 622-3304 | www.val-u-care.com


Val-U-Care
Fall Prevention & Anti-Wandering

Prevent Falls & Wandering

**CMS
 Compliant Fall
 Prevention**

Prevent Wandering With Anti-Wandering Door Systems

For Homes &
 Facilities of any
 size!



Affordable – No Wiring – No Electricians

MFR#	PRODUCT DESCRIPTION
TL-4005SYS	Single Door Bar Monitor System (1 year warranty)
TL-4012	Resident Wristband Transmitter (6 month warranty)
TL-4014W	Wall Mount Reset Button (1 year warranty)

Optional Additional Components:



TL-2025
 Magnetic Door Lock
 Will lock door when alarms



TL-4015
 Central Monitor

All-in-One

This "All in One" Anti-Wandering Door Bar System functions with just a door alert and a resident wristband. Simply mount the door bar by any exit or doorway in your care facility, plug it in, and place wristband on resident. When a resident wearing a wristband gets within 8 to 14 feet of the door bar or attempts to go through the doorway, the door bar will alert audibly and visually. An optional pager is available.

Simple Installation

1. Mount Door Bar
2. Activate Wristband Transmitter
3. Test



**CMS
Compliant for
Non-Restraint**

Prevent Falls & Wandering

Facility Grade Large Central Monitor Easy setup and works with large selection of wireless devices without any hard-wiring!

Large Facility Central Monitor
TL-4015



Caregiver Pager
TL-4016P



WIRELESS Monitor
TL-4016R

Soft Touch Button
NC-03



DepartAlert® Door Bar System
TL-4005SYS

Reset Button
4017-RB



Call Button
with Emergency Pull Cord
4017-CB



Emergency Call Alert
4017-ECA



Call Button
TL-4015NC



Corded
Bed Pad



Corded
Chair Pad



Corded
Floor Mat

Don't Let Them Walk Alone

Designed to wirelessly connect with up to 999 devices. Alerts caregiver when resident is in need of assistance and displays which device triggered the alarm. The Smart Station can identify alerting devices with a designated number. The Smart Station functions with a variety of components and accessories including, Anti-Wandering Door Systems, Emergency Call Button & Light Systems, LCD Pagers and Quiet Fall Prevention Alert to Pager Systems.

TL-4015 Key Features:

- Alert indicates which device triggered the alarm
- Sends alert to caregiver pager
- Receives signal from up to 999 wireless devices in the facility
- Batteries not included (4 x C)
- AC Adapter included (AC-04)
- Range: up to 300 feet



Call Us Today! (888) 622-3304 | www.val-u-care.com